

Bringing Vision to the Operating Room

The NUVO OR Suite – A Truly Integrated System


NUVO's integrated system of centrally-mounted surgical lights and monitors provides freedom of movement for optimal positioning and posture, plus maximum viewing flexibility during endoscopic procedures. With more than 2,000 lights in OR suites nationwide, our proven performance makes NUVO the light of choice in today's surgical environment.

Our patented, ceiling-mounted equipment features:

- Heavy-duty construction
- Drift-free positioning of lights and monitors
- Continuous rotation on three axes
- Multiple-image viewing
- Controls conveniently located on the sterile handles
- 136,500 lux of cool, shadow-free illumination from each light head



Multi-functional Sterile Handle



Easily Accessible Bulbs



Flat-folding Monitor Storage

NUVO Surgical Equipment – thoroughly-tested, industry-refined, proven results.

For more product information, call 800-663-1152 or visit www.nuvosurgical.com.



N U V O , I N C .

5368 Kuhl Road • Erie, Pennsylvania 16510-4703
(800) 663-1152 • (814) 899-4220 • FAX (814) 899-1410