



VERDÉ EX

LED EXAM LIGHT



BRINGING
VISION
TO THE OPERATING ROOM

PRODUCT
SPECIFICATIONS

MOUNTING OPTIONS



Mobile Stand



Ceiling Mount



Wall Mount

PRODUCT SPECIFICATIONS

Number of LEDs	19
Central illumination level at 6' (1m) distance	50,000 lux
Dim levels	0 / 50% / 75% / 100%
Power input (at dim level)	24 W (50%)
	32 W (75%)
	38 W (100%)
Color temperature	4,000 K
Color rendering index (R_a)	>94
Diameter of light field	6.5" (160 mm)
Depth of illumination	55" (1,400 mm)
Service life of light source	> 20,000 h
Power supply	100V 50-60 Hz
	120V 50-60 Hz
	220V 50-60 Hz
	230-240V 50-60 Hz
Size of light head	Ø 9.5" (240 mm)
Range of height adjustment of light head	52" (133 cm)
Protection category	IP 20



A member of the Medical Illumination International Group

5368 Kuhl Road • Erie, PA 16510-4703
 (800) 663-1152 • (814) 899-4220
 Fax: (814) 899-1410
www.nuvosurgical.com

NUVO reserves the right to make changes without notice in design, specifications, and models. The only warranty NUVO makes is the express written warranty extended on the sale or rental of its products.

© 2010 NUVO, Inc.
 ALL RIGHTS RESERVED.

