



VISTOR MS  
MINOR SURGERY LIGHTS



BRINGING  
VISION  
TO THE OPERATING ROOM

PRODUCT  
OVERVIEW




The applications of the VistOR MS Minor Surgery Light extend throughout the hospital into ambulatory surgery. This light combines unparalleled performance in a compact lightweight light head and delivers it at an extremely affordable price.

Whether your needs are for a surgical light in the OR or Trauma, a treatment light in the ER or simply a mobile exam and procedure light the VistOR MS can be configured to meet your needs.

The VistOR MS light is a true example of how high quality and high performance can be delivered in an extremely affordable manner. This light has created a new standard in the industry being available as Ceiling Mount, Wall Mount or Freestanding Mobile models thereby extending its range of applications beyond the hospital into ambulatory, small clinic and physician office settings.

NUVO ...  
sharing the vision



**16-inch Satellite Light Head**  
With light intensity of 72,000 lux, a 48-inch depth of field and a fully-adjustable light pattern, this light makes the perfect companion to any major surgical light. Its performance is further enhanced with the adjustable intensity settings and back-up bulb making it well suited to use for the most demanding procedures.



### 20-inch Surgical Light Head

The larger Light Head takes the performance to the next level, increasing the intensity to 102,000 lux while maintaining all the key attributes of the Satellite Light Head. The larger reflector provides even better shadow control making this light perfectly suited to virtually any procedure requiring surgical-quality lighting.



### V-cam Video Camera

The V-cam is a fully-integrated video camera system that delivers high-quality images with perfect clarity. The 10X Optical Zoom and 27X Digital Zoom, combined with its auto focus feature, allow the camera to be positioned optimally over the patient to gain the absolute best picture angle while maintaining the highest picture quality.

# CONFIGURATION OPTIONS



*20-inch Single Light Head shown*

## SINGLE LIGHT OPTION

The single light mount can be configured to have either the 20-inch or the 16-inch light head. Both lights feature continuous rotation of all the joints combined with ultra lightweight, easy-to-position counterbalance arms making these lights ideally suited to the needs of minor surgery.



*20-inch and 16-inch Light Head Combination shown*

## DUAL LIGHT OPTION

The dual light offers the flexibility of being set up with dual 16-inch or 20-inch light heads or a combination of both. The intensity of these lights combined with their excellent shadow control, variable pattern size and back-up bulb makes them an extremely popular choice for trauma and ambulatory surgery applications.



*Camera shown with a 20-inch and 16-inch Light Head Combination*

## DUAL LIGHT PLUS CAMERA OPTION

The flexible nature of the Vistor MS allows you to configure any combination of the single or dual light with the integrated camera system. The freedom of movement afforded by mounting the camera on a third arm provides optimum positioning during a case to capture still images for documentation or continuous video for teaching.



*A member of the  
Medical Illumination International Group*

5368 Kuhl Road • Erie, PA 16510-4703

(800) 663-1152 • (814) 899-4220

Fax: (814) 899-1410

[www.nuvosurgical.com](http://www.nuvosurgical.com)



S2Z00677

NUVO reserves the right to make changes without notice in design, specifications, and models. The only warranty NUVO makes is the express written warranty extended on the sale or rental of its products.

© 2007 NUVO ALL RIGHTS RESERVED.