



VISTOR PRO LED
EXAM AND PROCEDURE LIGHT






BRINGING
VISION
TO THE OPERATING ROOM



PRODUCT
SPECIFICATIONS

LED TECHNOLOGY IS NOW OUTSIDE OF THE OPERATING ROOM

CONFIGURATION OPTIONS



Available in single or dual ceiling mount, wall mount or mobile stand

The VISTOR PRO LED is another addition to the growing line of NUVO LED Lights. The state-of-the-art design provides unparalleled optics for ambulatory surgery centers, procedure rooms, trauma rooms, clinics and physician offices.

Control on/off and intensity settings from sterile handle.



The VISTOR PRO LED provides superior illumination and unmatched shadow control using 36 high efficiency ultracool LED's.

PRODUCT SPECIFICATIONS

Light Intensity	64,500 lux
Color Temperature	4500 K
Intensity Settings	5(50% – 100%)
Controls in Sterile Field	Yes
LED Life Efficiency	> 30,000 hrs.



*A member of the
Medical Illumination International Group*

5368 Kuhl Road • Erie, PA 16510-4703

(800) 663-1152 • (814) 899-4220

Fax: (814) 899-1410

www.nuvosurgical.com

